



Steps for Installing the Netripples Blood Bank System.

After downloading the software extract the downloaded software and find the below path of setup icon to run the set up.


Path : Open BBMS 5.0 setup folder, then open Netripples Blood Bank System, then open Product Configuration1 folder, open Release 1 folder, open Disk Images Folder, Open Netripples BBMS 5.0 Setup Folder, Find the Setup Icon and double click on it to run the setup.

Step 1: Click on the [Next >] button .On clicking [Next>] button it will open License Agreement screen. Read the EULA agreement carefully and click on the **I accept the terms in the license agreement** radio button .click the [Next>] button as shown below







Step 2: Enter User name and Organization and click on the Next Button




Step 3: Select Complete and Click on the Next Button to start the installation process.



Step 4: Click on the Install Button to start the installation process.



Step 5: Now the installation process may take few minutes while installing



Step 6: Click on the [Finish] button to complete the installation Process



After installation please follow the below step to run the application

Software Icons :

