





Steps for Installing the Netripples Dentos System.

After downloading the software extract the downloaded software and find the below path of setup icon to run the set up.


Path : Open BBMS 5.0 setup folder, then open Netripples Dentos System, then open Product Configuration1 folder, open Release 1 folder, open Disk Images Folder, Open Netripples Dentos 5.0 Setup Folder, Find the Setup Icon and double click on it to run the setup.

Step 1: Click on the [Next >] button .On clicking [Next>] button it will open License Agreement screen. Read the EULA agreement carefully and click on the **I accept the terms in the license agreement** radio button .click the [Next>] button as shown below







Step 2: Enter User Name and Organization name in the text fields and click on the [Next>] button.




Step 3: Select **Complete** radio button and click on the **[Next>]** button to start the installation process.



Step 4: Click on the **[Install]** Button to start the installation process.



Step 5: Now the installation process may take few minutes while installing



Step 6: Click on the [Finish] Button to complete the installation Process



After installation please follow the below step to run the application

Software Icons :

